



INTERNAL HEX Implant System



# infinite **OPPORTUNITIES IN IMPLANTOLOGY**

For almost 50 years ACE Surgical has been dedicated to dental surgical advancements. We continue to develop and manufacture the highest quality, state-of-the-art products at competitive prices; while keeping customer service at the core of our business.

Our latest development, the infinity Dental Implant System has been developed, engineered, and manufactured by some of the industry leading professionals in the USA. The infinity Dental Implant System allows you to place and restore our implants with the same confidence you get from your current system without the added expenses.

The infinity INTERNAL HEX Dental Implant System is committed to delivering a truly compatible dental implant solution to both you and your patients.

Infinity Implants—infinite innovation, endless opportunities in implantology.

4.7mm  
Implant

5.1mm  
Implant

3.7, 4.1, 4.7, 5.1mm I-HEX. DRILL SEQUENCE

FOR HARD BONE WITH COUNTERSINK & TAP

L (Drills & Tap)

Rev. A (DRILL & TAP SEQUENCE)

ACE Surgical Supply Co., Inc. • 1.800.441.3100 • [acesurgical.com](http://acesurgical.com)

1034 Pearl Street, Brockton, MA 02301 • E-Mail: [info@acesurgical.com](mailto:info@acesurgical.com)



Table of contents

Introduction ..... 3

Features and Benefits ..... 4

Implant Offerings ..... 6

Drilling Sequences ..... 10

Surgical Drills, Drivers, and Kit ..... 15

Prosthetics ..... 18

Surgical Protocols ..... 22

Prosthetic Protocols..... 26

Warranty ..... 30



**ACE Surgical Supply Co., Inc.**  
1034 Pearl Street, Brockton, MA 02301 - USA  
1-800-441-3100 • [www.acesurgical.com](http://www.acesurgical.com)

© Copyright 2015 ACE Surgical Supply Co., Inc.  
Reproduction in whole or part is specifically prohibited without written permission.



# Introduction

## **INTERNAL HEX IMPLANTS**

---

The logo for Infinity Dental Implant Systems is positioned in the lower half of the page. It features a large, light blue 'infinity' symbol (a horizontal figure-eight) above the word 'infinity' in a large, lowercase, sans-serif font. Below 'infinity' is the phrase 'Dental Implant Systems' in a smaller, uppercase, sans-serif font. The entire logo is set against a background of a large, light blue circle and a smaller, darker blue hexagon.

infinity

Dental Implant Systems

## infinity INTERNAL HEX Implants

The infinity INTERNAL HEX implant has been designed to work with some of the leading internal hex implants in the market today. The infinity INTERNAL HEX platform allows for a secure connecting interface between the implant and prosthetic. The precise manufacturing of all INTERNAL HEX components is what separates this system from other compatible dental implant systems.

The implants are manufactured in the USA from high-strength, bio-compatible titanium alloy (Ti 6Al-4V-ELI). The implant's roughened surface treatment allows for a greater surface area for osseointegration.

The engineered micro-grooves found on the collar of the implant are for maintaining crestal bone. The tapered thread design allows for immediate mechanical fixation, which in most instances, allows for early implant loading.

The available prosthetics options for the infinity INTERNAL HEX are designed to offer familiar options without the inflated costs. From a single tooth restoration to a full arch, infinity INTERNAL HEX prosthetics are easy to use, and restore functionality at a price that makes sense for your practice.



### INTERNAL HEX

Ø3.5mm Platform



Ø3.7mm  
Implant



Ø4.1mm  
Implant



### INTERNAL HEX

Ø4.5mm Platform



Ø4.7mm  
Implant




Ø5.1mm  
Implant

### COLOR-CODED

The entire INTERNAL HEX system is color-coded. From the drilling process, to the final restoration, we have color-coded all the components allowing you and your referrals to easily identify the implant platform and their correct components.

### RBM TEXTURED SURFACE



For over twenty years Resorbable Blast Media technology has allowed for improved osseointegration. RBM is a surface treatment designed to roughen the surface of an implant without leaving the residual embedded blast particles or debris in the treated substrate. To achieve the desired roughening, a biocompatible material with a suitable hardened particle size is used and then subsequently dissolved from the surface with a defined passivation treatment. The result is a rougher surface than the traditional acid etch-treatments, providing a greater surface area for osseointegration.

### INTERNAL HEX PLATFORM

The popular INTERNAL HEX allows for low-profile restorative options compatible with internal hex configured implants.

### MICRO-GROOVES

Placed at the neck of all INTERNAL HEX implants are micro-grooves that aid in maintaining crestal bone.

### TAPERED THREAD

Allows for immediate mechanical fixation.

### TRIPLE LEAD THREAD

For faster implant insertion.

### SHARP CUTTING FLUTES

Allows for easy implant insertion.

### LIFETIME WARRANTY

The ACE Surgical implant warranty program is designed to support all clinicians involved with the infinity Implant System.

(See page 30 for details)



#### FIXTURE-MOUNT/TRANSFER

Use as a fixture-mount or a transfer. The mount can be converted to a temporary abutment.

#### IMPLANT SUSPENSION CHAMBER

Color-coded for easy recognition. The suspension chamber securely cradles the implant and allows for easy removal once the cap is removed.

#### ALIGNMENT GUIDE

Implant is aligned with the flat, squared surface for correct buccal placement.

#### COVER SCREW

Included inside the vial's bottom cap.



#### STERILE PACKAGING

All implants come packed sterile with easy to read pull back seal with detailed labeling to expose the implant chamber.

#### CONTAINED IMPLANT

Vial packed implant suspension and easy to retrieve cover screw.

#### IMPLANT BOX

Implant symbol, size and color clearly distinguishes the INTERNAL HEX product. Side label includes implant specifications.

#### INSTRUCTIONS

Every implant complete with instructions for use enclosed.

#### PATIENT RECORD LABELS

Every implant comes complete with multiple patient labels for adhering to patient files.



**INTERNAL HEX****Ø3.7mm Tapered Implant**

Ø3.5mm Platform

Product No.	203708	203710	203711	203713	203716
Diameter	Ø3.7mm	Ø3.7mm	Ø3.7mm	Ø3.7mm	Ø3.7mm
Length	8mm	10mm	11.5mm	13mm	16mm

**IMPLANT  
COVER  
SCREW**Included with  
every Implant

20585001





## INTERNAL HEX

### Ø4.1mm Tapered Implant

#### Ø3.5mm Platform

Product No.	204106	204108	204110	204111	204113	204116
Diameter	Ø4.1mm	Ø4.1mm	Ø4.1mm	Ø4.1mm	Ø4.1mm	Ø4.1mm
Length	6mm	8mm	10mm	11.5mm	13mm	16mm



IMPLANT  
COVER  
SCREW

Included with  
every Implant



20585001

## INTERNAL HEX Healing Caps

#### Ø3.5mm Platform

- Uses a 1.25mm Hex Driver
- 30° Emergence Profile for Contoured Healing Caps
- Recommended Seating Torque 20 N-cm



Product No.	20135373	20135375	20135453	20135455	20135553	20135555
Platform	Ø3.5mm	Ø3.5mm	Ø3.5mm	Ø3.5mm	Ø3.5mm	Ø3.5mm
Collar Width	3.7mm	3.7mm	4.5mm	4.5mm	5.5mm	5.5mm
Collar Height	3mm Straight	5mm Straight	3mm Contoured	5mm Contoured	3mm Contoured	5mm Contoured

**INTERNAL HEX****Ø4.7mm Tapered Implant**

Ø4.5mm Platform

Product No.

204708

204710

204711

204713

204716

Diameter

Ø4.7mm

Ø4.7mm

Ø4.7mm

Ø4.7mm

Ø4.7mm

Length

8mm

10mm

11.5mm

13mm

16mm

**IMPLANT  
COVER  
SCREW**Included with  
every Implant


20585501

## INTERNAL HEX

### Ø5.1mm Tapered Implant

Ø4.5mm Platform

Product No.	205106	205108	205110	205111	205113	205116
Diameter	Ø5.1mm	Ø5.1mm	Ø5.1mm	Ø5.1mm	Ø5.1mm	Ø5.1mm
Length	6mm	8mm	10mm	11.5mm	13mm	16mm



IMPLANT  
COVER  
SCREW

Included with  
every Implant



20585501

## INTERNAL HEX Healing Caps

Ø4.5mm Platform


- Uses a 1.25mm Hex Driver
- 30° Emergence Profile for Contoured Healing Caps
- Recommended Seating Torque 20 N-cm




Product No.	20145473	20145475	20145553	20145555	20145653	20145655
Platform	Ø4.5mm	Ø4.5mm	Ø4.5mm	Ø4.5mm	Ø4.5mm	Ø4.5mm
Collar Width	4.7mm	4.7mm	5.5mm	5.5mm	6.5mm	6.5mm
Collar Height	3mm Straight	5mm Straight	3mm Contoured	5mm Contoured	3mm Contoured	5mm Contoured

# infinity **INTERNAL HEX 3.7mm IMPLANT** Drilling Sequences

## Soft Bone Drilling Sequence\*




## Dense Bone Drilling Sequence\*



\* Note: Instruments shown are drilling for a 11.5mm implant. Adjust accordingly for other implant depths.

**Caution:** The surgical twist drills are approximately 1mm longer in length than what is described.


Twist Drills and  
Depth Gauge  
marking Guide




# infinity **INTERNAL HEX 4.1mm IMPLANT**

## Drilling Sequences

### Soft Bone Drilling Sequence\*




### Dense Bone Drilling Sequence\*



\* Note: Instruments shown are drilling for a 11.5mm implant. Adjust accordingly for other implant depths.  
**Caution:** The surgical twist drills are approximately 1mm longer in length than what is described.


Twist Drills and  
Depth Gauge  
marking Guide




# infinity **INTERNAL HEX** 4.7mm IMPLANT

## Drilling Sequences

### Soft Bone Drilling Sequence\*




### Dense Bone Drilling Sequence\*



\* Note: Instruments shown are drilling for a 11.5mm implant. Adjust accordingly for other implant depths.  
**Caution:** The surgical twist drills are approximately 1mm longer in length than what is described.

Twist Drills and  
Depth Gauge  
marking Guide







# infinity **INTERNAL HEX 5.1mm IMPLANT**

## Drilling Sequences


### Soft Bone Drilling Sequence\*



### Dense Bone Drilling Sequence\*



Twist Drills and  
Depth Gauge  
marking Guide



\* Note: Instruments shown are drilling for a 11.5mm implant. Adjust accordingly for other implant depths.

**Caution:** The surgical twist drills are approximately 1mm longer in length than what is described.



# infinity **INTERNAL HEX** Surgical Drills, Drivers, and Kit



The image displays a variety of surgical instruments from the ACE Infinity Internal Hex system. In the foreground, a grey torque wrench is visible, labeled "Implant Torque Wrench". Behind it, several surgical drills and drivers are arranged on a blue surface. The drills are color-coded: blue for "Regular Drill", green for "Tap", and pink for "Countersink". The drivers are also color-coded: blue for "Drivers", green for "Countersink", and pink for "Countersink". The background is a solid blue color with a large, faint hexagonal graphic.

### CERTIFIED PRECISION

Manufactured under strictly controlled machining and certified to ISO 13485 regulations.

### BUILT TO LAST

Manufactured from hardened, high-strength, corrosive resistant surgical drilling materials. Engineered for multiple use.

### PROVEN DESIGN


Proven twist drill geometry allows efficient bone chip removal.



### STEP DESIGN

The steps on the larger drills fit the diameter of the previously drilled osteotomy to aid in alignment.

### LASER MARKED

Precise laser markings allows for accurate drill identification and implant depth drilling.





\* NEW Dense Bone Drill for 5.1 implant added. The surgical tray will be updated to accommodate this drill in early 2016



Ø 4.7mm Bone Tap  
04-9523-01

Ø 5.1mm Bone Tap  
04-9525-01



Ø 3.7mm Countersink  
04-9535-01




Ø 4.1mm Countersink  
04-9535-04



Ø 4.7mm Countersink  
04-9535-02



Ø 5.1mm Countersink  
04-9535-03



Complete infinity  
Internal HEX Surgical Kit  
20-9910-00



Ø 2.8mm Depth Gauge  
20-9920-01



Ø 2.0/2.8mm Paralleling Pins  
20-9730-02



2.5mm Hex CA Driver Tool  
04-9770-01



2.5mm Hex Ratchet Driver Tool  
20-1560-01



1.25mm Hex CA Driver  
04-1510-02



1.25mm Hex Ratchet Driver  
20-1510-02



1.25mm Hex Hand Driver  
20-1575-02



Drill Extender  
04-9599-01



22mm Square Ratchet Driver Extension  
20-1557-11



Open-End Offset Wrench  
20-1610-01



Implant Torque Ratchet - 20-45 N-cm  
003-5245



# infinity

## INTERNAL HEX Prosthetics



INTERNAL HEX Ø3.5  
PLATFORM PROSTHETICS  
Fits Ø3.7mm and Ø4.1mm Implants

### Implant Analog



Product No.

20735

### Closed and Open Tray Transfers

Impression Copings w/Retaining Screw



- Implant Level Transfer
- Includes Retaining Screw
- Uses a 1.25mm Hex Driver
- Recommended Seating Finger Pressure

Product No.

2063545

2053545

Collar Width

4.5mm

4.5mm

Tray Type

Closed Tray

Open Tray

### Prepable Abutments

with Retaining Screw



- Includes Retaining Screw
- Uses a 1.25mm Hex Driver
- Recommended Seating Torque 30 N-cm
- 1mm Prepable Score Lines

Product No.

20235452

20235454

Collar Width

4.5mm

4.5mm

Collar Height


2mm

4mm

Total Height

7mm

9mm



CW = Collar Width  
CH = Collar Height  
TH = Total Height

### UCLA



- Includes Retaining Screw
- Uses a 1.25mm Hex Driver
- Recommended Seating Torque 30 N-cm

Product No.

2093545

2083545

Collar Width

4.5mm

4.5mm

Engagement

Non-Engaged

Engaged





INTERNAL HEX Ø4.5mm  
PLATFORM PROSTHETICS  
Fits Ø4.7mm and Ø5.1mm Implants

### Implant Analog



Product No.

20745

### Closed and Open Tray Transfers

Impression Copings w/Retaining Screw



- Implant Level Transfer
- Includes Retaining Screw
- Uses a 1.25mm Hex Driver
- Recommended Seating Finger Pressure

Product No.

2064555

2054555

Collar Width

5.5mm

5.5mm

Tray Type

Closed Tray

Open Tray

### Prepable Abutments

with Retaining Screw



- Includes Retaining Screw
- Uses a 1.25mm Hex Driver
- Recommended Seating Torque 30 N-cm
- 1mm Prepable Score Lines

Product No.

20245552

20245554

Collar Width

5.5mm

5.5mm

Collar Height


2mm

4mm

Total Height

7mm

9mm



CW = Collar Width  
CH = Collar Height  
TH = Total Height

### UCLA



- Includes Retaining Screw
- Uses a 1.25mm Hex Driver
- Recommended Seating Torque 30 N-cm

Product No.

2094555

2084555

Collar Width

5.5mm

5.5mm


Engagement

Non-Engaged

Engaged



INTERNAL HEX Ø3.5 Platform  
**ZEST LOCATOR® ABUTMENT**  
**Fits Ø3.7mm and Ø4.1mm Implants**




Product No.	206-8662	206-8663	206-8664	206-8665	206-8625
Diameter	Ø3.5mm	Ø3.5mm	Ø3.5mm	Ø3.5mm	Ø3.5mm
Collar Height	1.0mm	2.5mm	3.5mm	4.5mm	5.5mm

- Recommended Seating Torque 30 N-cm



INTERNAL HEX Ø4.5mm Platform  
**ZEST LOCATOR® ABUTMENT**  
**Fits Ø4.7mm and Ø5.1mm Implants**



Product No.	206-8672	206-8673	206-8674	206-8675	206-8627
Diameter	Ø4.5mm	Ø4.5mm	Ø4.5mm	Ø4.5mm	Ø4.5mm
Collar Height	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm

- Recommended Seating Torque 30 N-cm

## ZEST LOCATOR® COMPONENTS

### Fits All INTERNAL HEX Implants

#### REPLACEMENT MALE - SETS

Each Set includes:

- (1) Denture Cap with Black Processing Male
- (1) White Block-Out Spacer
- (3) Nylon Retention Males: 1 each of the Dual or Extended Range Nylon Retention Males

#### Dual Retention Replacement Male

Two retention surfaces, the outer ring and the center nipple allows for maximum hold. Good for use with divergent implants up to 10° (accommodates up to 20° between 2 or more implants).

#### Extended Range Replacement Male

Single retention surface using only the outer ring with the inner nipple removed. Good for use with divergent implants up to 20° (accommodates up to 40° between 2 or more implants).

#### REPLACEMENT PROCESSING PACKS

- Available in various retention strengths
- Sold in packages of 4



Zest LOCATOR® Replacement  
Dual Retention Male Sets

2 sets / pkg - 206-85192  
10 sets / pkg - 206-851910



Zest LOCATOR® Replacement  
Extended Range  
Male Sets

2 sets / pkg - 206-85402  
10 sets / pkg - 206-854010




Color	CLEAR	PINK	BLUE	GREEN	ORANGE	RED
Product No.	206-8524	206-8527	206-8529	206-8547	206-8915	206-8548
Retention	Standard	Light	Extra Light	Standard	Light	Extra Light
Lbs.	5 lbs.	3 lbs.	1.5 lbs.	3-4 lbs.	2 lbs.	.5-1.5 lbs.
Grams	2268	1361	680	1361-1814	907	226-680
	Male Processing Pack (Dual Retention)			Male Processing Pack (Extended Range)		

#### CHAIRSIDE DENTURE PREP & POLISH KIT

### For Installation of Locator® Denture Caps

- A comprehensive kit consisting of Recess, Trim, Undercut, Grind and Vent Burs, as well as a Polisher, all of which are designed to address the most frequent overdenture preparation requirements.
- A unique to the market CHAIRSIDE Recess Bur, specifically designed by ZEST Anchors, that quickly and easily prepares the exact size recess for the LOCATOR and SATURNO Denture Attachment Housings (Denture Caps)



Zest CHAIRSIDE  
Denture Prep  
& Polish Kit  
206-09582



Zest LOCATOR®  
Torque Wrench  
Contra Angle Driver  
206-8914

Zest LOCATOR® Male  
Installation / Removal Tool  
This end for LOCATOR® Abutments

#### ZEST LOCATOR® LAB COMPONENTS



Zest LOCATOR®  
White Block-Out  
Spacer - 20/pkg  
206-8514



Zest LOCATOR®  
Transfer - 4/pkg  
206-8505



Zest LOCATOR®  
Analog - 4/pkg  
206-8516



Zest LOCATOR® Male  
Installation / Removal Tool  
206-8393



# infinity INTERNAL HEX Surgical Protocol




## Implant Site Preparation

### INTERNAL HEX Implant Placement Protocol

Drilling depth shown for 11.5mm(L) implant.

Twist Drills and  
Depth Gauge  
marking Guide



#### Step 1 - Determine Tissue Depth

- Prior to drilling, use a tissue probe to determine the soft tissue depth.

#### Step 2 & 3 - Preparing the Surgical Site

- Make a mesio-distal incision along the buccal side of the alveolar crest through the mucoperiosteum and attached gingiva to the bone and reflect the flap.

#### Step 4 - Marking the Implant Site

- Using a handpiece at 800 to 1200rpm and copious irrigation, mark the initial osteotomy using the Ø1.6mm round drill.

#### Step 5 - Drilling the Pilot Hole

- Drill the initial pilot osteotomy using the Ø2.0mm twist drill.
- Drill to appropriate implant length.

#### Step 6 - Check Parallel Alignment

- Confirm the appropriate drill angle with the 2mm paralleling pin.
- Osteotomy corrections can be made during the following drilling step.

#### Step 7A - Initial Implant Drill - Ø2.8mm/Ø2.3mm

- Select the Ø2.8mm twist drill with the Ø2.3mm step and drill to the appropriate length.


#### Step 7B - Measure Drill Depth

- Check to make sure you have achieved the appropriate depth by placing the depth gauge into the osteotomy.

**Ø3.7 NOTE - If placing a Ø3.7mm Implant:**

**FOR SOFT BONE**, this is the final step. Proceed to "Loading Implant Onto Driver" on page 24.

**FOR DENSE BONE**, proceed to "Additional Steps for Dense Bone" on page 23.



Continued on following page



## Implant Site Preparation - Cont.

### INTERNAL HEX Implant Placement Protocol

Drilling depth shown for 11.5mm(L) implant.

#### Step 7C - Twist Drill - Ø3.4mm/Ø2.8mm

- Select the Ø3.4mm twist drill with the Ø2.8mm step and drill to the appropriate length.

#### Ø4.1 NOTE - If placing a Ø4.1mm Implant:

**FOR SOFT BONE**, this is the final step. Proceed to "Loading Implant Onto Driver" on page 24.

**FOR DENSE BONE**, proceed to "Additional Steps for Dense Bone" below, on this page.

#### Step 7D - Twist Drill - Ø3.8mm/Ø3.4mm

- Select the Ø3.8mm twist drill with the Ø3.4mm step and drill to the appropriate length.

#### Ø4.7 NOTE - If placing a Ø4.7mm Implant:

**FOR SOFT BONE**, this is the final step. Proceed to "Loading Implant Onto Driver" on page 24.

**FOR DENSE BONE**, proceed to "Additional Steps for Dense Bone" below, on this page.




#### Step 7E - Twist Drill - Ø4.4mm/Ø3.8mm

- Select the Ø3.8mm twist drill with the Ø3.4mm step and drill to the appropriate length.

#### Ø5.1 NOTE - If placing a Ø5.1mm Implant:

**FOR SOFT BONE**, this is the final step. Proceed to "Loading Implant Onto Driver" on page 24.

**FOR DENSE BONE**, proceed to "Additional Steps for Dense Bone" below, on this page.



## Additional Steps for Dense Bone

### INTERNAL HEX Implant Placement Protocol

Drilling depth shown for 11.5mm(L) implant.

#### Step 8- Twist Drill - Use the Next Larger Size:

- Ø3.7mm implant, use the Ø3.4mm/2.8mm step, Ø4.1mm implant, use the Ø3.8mm/3.4mm step, Ø4.7mm implant, use the Ø4.4mm/3.8mm step, Ø5.1mm implant, use the Ø4.8mm/4.4mm step,
- Drill to the appropriate length.



#### Step 9- Countersink - Match the Implant Diameter:

- Select the matching diameter countersink and drill to the top of the laser marking.

#### Step 10- Bone Tap - match the Implant Diameter:

- Select the matching diameter bone tap and using slow speed (20rpm) tap the osteotomy to the bottom of the laser marking as shown.

Please refer to pages 10 - 13 for the complete drilling sequences.



Instruments shown are for reference only. Follow Dense Bone Additional Steps 8, 9, and 10 for actual instrument sizes.




## Loading Implant Onto Driver

### INTERNAL HEX Implant Placement Protocol

#### Step 1 - Implant Packaging

- Remove the Tyvek® sealed blister package from the outside box. Open this blister package over a sterile field by locating and pulling off the sealed Tyvek® lid. With the blister package open, drop the sterile inner implant vial onto the sterile field.




#### Step 2 - Break Sterile Strip - Remove Cap

- Open the sterile implant vial by breaking the tamper proof tape and removing the cap using a gentle twisting and pulling motion.

#### Step 3 - Loading the Driver

- Using either the 2.5mm hex ratchet or contra-angle driver tool, carefully rotate and engage the hexagon of the implant and driver tool in a downward direction.
- Pick up the implant body from the vial.
- Use caution when bringing implant to the surgical site.



## Placing Implant Into Osteotomy

### INTERNAL HEX Implant Placement Protocol

#### Step 1 - Placement Setup and Speed

- With the implant attached to the driver tool, insert the implant into the prepared osteotomy using an implant placement speed of 20 rpm and a maximum insertion torque of 45 N-cm.
- If the implant does not seat at 45 N-cm, back the implant out of the osteotomy, place the implant back into the sterile vial, then refer to the dense bone protocols (Step 7B-E on page 22-23).

\* External irrigation may be used to minimize heating during this process.

#### Step 2 - Seating the Implant





- Drive the implant until fully seated and flush with the crest of the bone ridge.
- Do not exceed 45 N-cm.

#### Step 3 - Disengage Driver

- After placement, disengage the hand or contra-angle implant driver tool from the implant and fixture-mount by pulling off the tool in a straight upward direction.

#### Step 4 - Remove Fixture-Mount

- Using the offset wrench to keep the fixture-mount from turning, remove the retaining screw from the top of the fixture-mount using the 1.25mm hex hand driver. This will remove the fixture-mount from the implant.








## Two-Stage Healing Option

### INTERNAL HEX Implant Placement Protocol



#### Step 1 - Loading the Implant Cover Screw

- The implant cover screw is located inside the bottom cap of the implant vial. Using one of the 1.25mm hex drivers, turn counter-clockwise to release the screw from the cap.



#### Step 2 - Placing the Implant Cover Screw

- Carefully place the implant cover screw into the implant and securely tighten into place.



#### Step 3 - Closing Tissue

- Suture the tissue over the cover screw using standard surgical protocols.



## Single-Stage Healing Option

### INTERNAL HEX Implant Placement Protocol

#### Step 1 - Healing Cap

- Using one of the 1.25mm hex drivers, place the appropriate healing cap into the implant using the recommended seating torque of 20 N-cm and secure tightly into place.



#### Step 2 - Closing Tissue

- Suture the tissue around the healing cap using standard surgical protocols.





## Single-Stage Temporary Preable Abutment Healing Option

### INTERNAL HEX Implant Placement Protocol

#### Step 1 - Remove Fixture-Mount

- Using the offset wrench to keep the fixture-mount from turning, remove the retaining screw from the top of the fixture-mount using the 1.25mm hex hand driver. This will remove the fixture-mount from the implant.





#### Step 2 - Attach Fixture-Mount to Analog

- Attach fixture-mount to its corresponding color-coded implant analog.


#### Step 3 - Remove Fixture-Mount

- Using a circular saw remove the top colored section from fixture-mount through groove, creating a straight abutment.
- As required, prepare the temporary abutment to the appropriate occlusal height.



#### Step 4 - Attach the Temporary Abutment

- Use a 1.25 driver to attach abutment to the implant and torque the temporary abutment to 30 N-cm.



#### Step 5 - Closing Tissue

- Suture the tissue around the temporary abutment using standard surgical protocols.

#### Step 6 - Placing the Temporary Crown

- After prepping closing the tissue around the temporary abutment, cement the temporary crown to the abutment.

# infinity INTERNAL HEX Prosthetic Protocol



## INTERNAL HEX Implant Prosthetic Protocol Implant Level Impression Technique - Closed Tray

### Step 1 - Removing the Healing Cap or Cover Screw

- Using one of the 1.25mm hex drivers, remove the healing cap or cover screw from the implant.

Note: In the case of bony overgrowth, use a curette to carefully clear the surface of the implant to provide direct access for transfer seating.

### Step 2 - Place the Closed Tray Impression Transfer

- Use finger pressure to seat the closed tray transfer with transfer retaining screw using one of the 1.25mm hex drivers.
- A radiograph may be taken to verify proper seating of the closed tray impression transfer.

### Step 3 - Impression Preparation

- Use cotton pellet or material to protect the head of the retaining screw. Fill the access hole to prevent impression material from entering.
- Using impression material, encase the transfer.

### Step 4 - Take Impression


- Place the loaded impression tray into the mouth, apply pressure and allow the impression materials to set.

### Step 5 - Removing Transfer

- Once impression tray is lifted, unblock the transfer retaining screw opening. Using the 1.25mm hex driver, remove the transfer.
- Remove the impression transfer and materials protecting the head of the retaining screw.
- Thread the impression transfer onto the corresponding implant analog.
- Send the complete implant analog assembly and the impression to dental laboratory for model fabrication.

### Step 6 - Replace Healing Cap

- Re-seat the healing cap with a 1.25mm hex driver using the recommended seating torque of 20 N-cm.





## INTERNAL HEX Implant Prosthetic Protocol

### Implant Level Impression Technique - Open Tray

#### Step 1 - Removing the Healing Cap or Cover Screw

- Using one of the 1.25mm hex drivers, remove the healing cap or cover screw from the implant.

Note: In the case of bony overgrowth, use a curette to carefully clear the surface of the implant to provide direct access for transfer seating.

#### Step 2 - Placing the Open Tray Impression Transfer

- Seat the open tray transfer with transfer retaining screw using the 1.25mm hex driver.

#### Step 3 - Prepare An Opening in Tray

- Mark and then cut a Ø5-10mm opening in the top of the impression tray at the exact area where the transfer will protrude.

#### Step 4 - Impression Preparation

- Using impression material, encase the transfer.

#### Step 5 - Take Impression


- Place the loaded impression tray into the mouth, apply pressure, and allow both impression materials to set.

#### Step 6 - Removing Tray and Transfer

- Using one of the 1.25mm hex driver, carefully remove just the retaining screw through the access hole, leaving the open tray transfer in the impression.
- Once impression has set, passively remove the impression with the remaining impression transfer embedded.
- Send the impression along with the laboratory components to the dental laboratory for model fabrication.

#### Step 7 - Replace Healing Cap

- Clean the implant area and healing cap.
- Re-seat the healing cap with a 1.25mm hex driver using the recommended seating torque of 20 N-cm.







## INTERNAL HEX Implant Prosthetic Protocol

### Prepable Abutment Impression Technique

#### Step 1 - Removing the Healing Cap or Cover Screw

- Using one of the 1.25mm hex drivers, remove the healing cap or cover screw from the implant.

Note: In the case of bony overgrowth, use a curette to carefully clear the surface of the implant to provide direct access for final abutment seating.

#### Step 2 - Preparing the Prepable Abutment

- Seat the abutment with the included retaining screw and mark for preparations. Mark the abutment's final height and gingival contour.

#### Step 3 - Abutment Preparation

- Remove abutment and attach to implant analog and make any necessary adjustments.
- Re-seat the abutment with the included retaining screw and torque to 30 N-cm.

#### Step 4 - Impression Preparation






- Use cotton pellet or material to protect the head of the retaining screw. Fill the access hole to prevent impression material from entering the abutment.
- Using impression material, encase the abutment.

#### Step 5 - Take Impression

- Place the loaded impression tray into the mouth, apply pressure, and allow impression material to set.
- Send impression to laboratory.

#### Step 6A, 6B - Temporary Process

- A temporary crown can be made and cemented to the prepped abutment or left in place while the final prosthetic is being made.





## INTERNAL HEX Implant Prosthetic Protocol **LOCATOR® Abutment**

### Step 1 - Removing the Healing Cap or Cover Screw

- Using one of the 1.25mm hex drivers, remove the healing cap or cover screw from the implant.
- Note: In the case of bony overgrowth, use a curette to carefully clear the surface of the implant to provide direct access for transfer seating.

### Step 2 - Placing the LOCATOR® Abutment

- Seat the LOCATOR abutment by hand using finger pressure.

### Step 3A, 3B - Torque the LOCATOR® Abutment


- Once seated on the implant, torque the LOCATOR Abutment to 30 N-cm using the brass end of the LOCATOR male installation tool with a standard 1.25mm hex driver.

### Step 4 - Impression Preparation

- Place the LOCATOR Impression Transfers onto the abutments.

### Step 5A, 5B, 5C - Take Impressions

- Using impression material, cover the transfers.
- Place the loaded impression tray into the mouth, apply pressure, and allow impression material to set.
- Remove impression and send along with laboratory components to the laboratory.



# infinity Implant Warranty

The ACE Surgical implant warranty is designed to support our customers that place infinity Dental Implant Systems. This warranty will address all aspects of the procedure including the surgical and restorative components.

To make a warranty claim:

- Call our Customer Experience department to report complaint and obtain an implant warranty form.
- Identical or similar implant/product will be sent and invoiced.
- Return the contaminated implant packaged sterile with the completed implant warranty form.
- Upon receipt of the implant, ACE will credit the returned implant.

The warranty will cover:

- A Same day spinners**
- B Post-placement, pre-loading implant failure**
- C Post-loading failure**
- D Post-loading failure using ACE restorative components on a different implant manufactured implant**

## **A Same day spinners**

ACE will replace any implant that for any reason could not be placed at the time of surgery.

## **B Post-placement, pre-loading implant failure\***

ACE will provide the following:

- An identical implant.
- A \$150 product credit on any ACE regenerative products, sutures or instrumentation needed to graft the site.

## **C Post-loading failure**

ACE will replace any implant that fails post-loading.

ACE will provide the following\*:

- ACE restorative components.
- A \$150 product credit on any ACE regenerative products, sutures or instrumentation needed to graft the site.
- One additional ACE implant and corresponding restorative components.

## **D Post-loading failure using ACE restorative components on a different implant manufactured implant\***

ACE will provide the following:

- Replacement of the restorative components.
- A similar implant as to the one used in the surgery.

\*Submit a radiograph showing the compromised site.



Please contact your ACE representative with any warranty related issues.

**ACE Surgical Supply Co., Inc.**

1034 Pearl Street, Brockton, MA 02301 - USA

Tel. 1-800-441-3100 Fax 1-800-583-3150

[www.acesurgical.com](http://www.acesurgical.com)





**ACE SURGICAL SUPPLY CO., INC.**

1034 Pearl Street, Brockton, MA 02301 - USA

1-800-441-3100 • [www.acesurgical.com](http://www.acesurgical.com)